



THUASNE

KITBIFLEX - KITBIFLEX Hôpital

| | | |
|----|--|----|
| fr | Kit de 2 bandes lavables de compression élastique, à allongement court..... | 4 |
| en | Kit of 2 short-stretch washable elastic compression bandages | 6 |
| de | Kit mit 2 waschbaren, elastischen Kurzzugkompressionsbinden | 8 |
| nl | Set van 2 wasbare zwachtels met korte rek | 10 |
| it | Kit di 2 bende elastiche di compressione lavabili, ad allungamento limitato..... | 12 |
| es | Equipo de 2 vendas lavables de compresión elástica, de prolongación corta..... | 14 |
| pt | Kit de 2 faixas laváveis de compressão elástica, de extensão curta | 16 |
| da | Sæt med 2 vaskbare elastiske kompressionsbind, med kort stræk | 18 |
| fi | Pakkaus, joka sisältää 2 pestävää vähäelastista venyvää puristussidettä | 20 |
| sv | Kit om 2 tvättbara, elastiska, kortsträckta kompressionsbandage..... | 22 |
| el | Kit 2 πλενόμενων ελαστικών επιδέσμων συμπίεσης, χαμηλής ελαστικότητας..... | 24 |
| cs | Sada dvou omyvatelných kompresivních elastických obinadel, krátkotažná obinadla..... | 27 |
| pl | Zestaw 2 elastycznych bandaży o krótkim naciągu, nadających się do prania..... | 29 |
| lv | Divu mazgājamu elastīgo kompresijas saišu komplekts, ar īsu izstiepumu | 31 |
| lt | Dviejų plaunamų, tamprių kompresinių, nedaug išsitempiančių juostų rinkinys | 33 |
| et | Komplekt 2 pestava elastse survesidemega, väheveniv..... | 35 |
| sl | Komplet dveh pralnih elastičnih kompresijskih povojev z majhnim raztezanjem | 37 |
| sk | Súprava 2 elastických kompresívnych krátkotažných obvazov, ktoré možno prať | 39 |
| hu | 2 darab, mosható, rövid megnyúlású, rugalmas kompressziós pólya..... | 41 |
| bg | Комплект от 2 еластични стягащи бинтове, годни за пране, слабо разтегливи | 43 |
| ro | Kit de 2 benzi lavabile de compresie elastică, cu întindere scurtă..... | 45 |
| ru | Набор из 2 эластичных компрессионных бинтов малой растяжимости. Изделия пригодны для стирки. | 47 |
| hr | Komplet dva periva elastično-kompresivna zavoja s kratkim rastezanjem | 50 |
| zh | KIT內含2條可清洗低延展性彈性繃帶..... | 52 |
| ja | 洗浄可能弾性圧迫包帯2個からなるキット、ショートストレッチ..... | 53 |
| ko | 세탁 가능한 저탄력 압박 붕대 2개로 구성된 키트..... | 55 |
| ar | طاقم من رباطين ضاغطين قابلين للغسل مرتين، لهما استطالة قصيرة..... | 59 |

I

II

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(1) Haute Autorité de Santé (France)

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en

KIT OF 2 SHORT-STRETCH WASHABLE ELASTIC COMPRESSION BANDAGES

Description

The KITBIFLEX consists of two bandages:

Bandage No. 1 (I): Compression bandage which acts to provide support and protect the skin - short-stretch elastic bandage (stretch of 10 to 100%) - (polyester microfibre + viscose + elastane + polyamide).

Bandage No. 1 protects the skin and bony areas and supports the limb.

Bandage No. 2 (II): Compression bandage - short-stretch elastic bandage (maximum stretch of 10 to 100%) - (viscose - polyester - polyamide - elastane). Bandage No. 2 provides the required pressure level.

The KITBIFLEX is available in 3 sizes depending on ankle girth:

| Kit size | Ankle size | Bandage 1 length | Bandage 2 length |
|----------|------------|------------------|------------------|
| Size 1 | < 25 cm | 4 m | 3.5 m |
| Size 2 | 25-30 cm | 4.5 m | 4 m |
| Size 3 | > 30 cm | 5 m | 4.5 m |

For treatment of venous ulcers, refer to the ankle girth.

For application on the upper limb, the choice of size is left to the judgement of the health care professional depending on the desired bandage length.

Properties/Mechanism of action

The KITBIFLEX provides compression, an essential component in the healing process of venous ulcers and in the treatment of chronic oedema or lymphoedema.

The bandages are washable and reusable.

For venous ulcers:

The combined action of the 2 KITBIFLEX bandages guarantees the required pressure level recommended by the HAS (French Health Authority) (> 36 mmHg).

These bandages can be left in place on the limb day and night, for a maximum of 7 days.

At the start of treatment of an ulcer, the bandage must be replaced and/or repositioned each time that local dressings are changed, i.e. as often as the exudate level requires, or whenever the volume of the limb changes (oedema).

It is recommended that treatment be continued until the venous ulcer has healed completely.

The continued use of a compression device (stockings or bandages) after healing is recommended to prevent recurrences.

For the treatment of chronic oedema:

Short-stretch dry bandages can be used as an alternative to compression stockings.

For the treatment of lymphoedema:

The KITBIFLEX short-stretch dry bandages deliver a therapeutic pressure level of over 36 mmHg as recommended for the treatment of lymphoedema. These bandages can be left in place day and night, according to the recommendations of the health care professional.

During treatment of lymphoedema, the bandage must be replaced and/or repositioned each time the volume of the affected limb changes.

The duration of treatment is recommended by the health care professional according to disease progression.

Use

The KITBIFLEX is intended for use on the upper and lower limbs of adult or paediatric patients suffering from at least one of the indications described below.

Indications

Venous ulcer of the lower limb (CEAP classification stage C6).

Chronic venous oedema.

Lymphoedema.

Instructions for use/Application

For the lower limb, it is preferable to apply the KITBIFLEX bandages in the morning after the patient gets up or after elevating the patient's legs for an hour.

Before applying the bandages:

Examine the shape of the limb; if there are any deformed or bony areas, use the protective foam bandages or Varico pads to make the shape of the limb uniform.

For application below the knee: apply the bandages in the following order keeping the foot at a right angle. To anchor the bandage on the foot, apply just enough pressure to hold the bandage in place. **Only use the markings from the ankle upwards.**

Bandage No. 1 (see application diagram I)

① - Position the end of the bandage flat against the base of the toes (softest side against the skin and calibrated side visible).

② - Fully cover the heel.

③ - Securing from the bottom.

④ - Securing from the top.

⑤ - Apply using a spiral method completely covering the square. Finish the bandage 2 fingers below the back of the knee.

Bandage No. 2 (see application diagram II)

① - Position the end of the bandage flat against the base of the toes allowing bandage No. 1 to show slightly (side with markings visible).

Any excess bandage should not be wrapped around the limb again since this could create a pressure peak. Either cut off the excess bandage or allow it to hang loosely lengthwise on the bandage.

To roll up the bandage correctly after use, follow the direction shown on the diagram (squares inside).

The KITBIFLEX (does not apply to the KITBIFLEX Hôpital):

How to use the bandage fasteners correctly (see Figures A and B):

Ⓐ - To apply the short fasteners:

• Place one side of the fastener on the end of the bandage to be secured in place.

• Pull gently on the fastener after application so as to apply slight tension to the bandage before securing the other side of the fastener in place. The attachment will be all the more effective if the bandage is slightly stretched.

• Press gently on the entire fastener to anchor it firmly on the bandage.

To use the long fastener:

If necessary, to limit the risk of slippage, wrap the long fastener around the top of the calf (see diagram Ⓒ), with no tension. The fastener must pass right round the limb and attach to itself with the Biflex logo visible.

Ⓑ - To remove the fasteners:

• Lift slowly by one edge until the fastener detaches (as if removing a dressing). It is possible that a few fibres or fluff could be left after the fastener is detached, but this does not affect the quality of the bandage or reduce its effectiveness.

5 short fasteners and 1 long fastener are supplied (does not apply to the KITBIFLEX Hôpital) in order to allow for replacement if they are lost.

Markings:

The markings help control the level of compression of the bandage.

The bandage is correctly stretched when each rectangle is transformed into a square.

The bandage has been applied correctly when the square is completely covered (following the blue line along the bandage).

The marking system and the application technique are identical for both bandages, thereby guaranteeing reproducibility of application and compliance with the appropriate pressure level required for treatment.

All bandages can be repositioned to adapt to changes in limb volume.

Contraindications

• Arterial pathologies: peripheral arterial disease (PAD), arterial ulcer or mixed ulcer with a predominantly arterial component (ankle-brachial index: ABI < 0.6 (according to the French Health Authority) or > 1.3).

• Patients suffering from advanced diabetic microangiopathy phlegmasia cerulea dolens with arterial compression, or septic thrombosis.

• Extra-anatomic bypass.

• Uncontrolled heart failure.

• Pre-existing or infected dermatological conditions (abscesses, boils, infected dermatitis).

• Known allergy to one of the components.

• Advanced peripheral neuropathy.

• Uncontrolled diabetes.

Specifically when worn on an upper limb:

• Brachial plexus pathology.

• Vasculitis of the extremities.

Precautions

• A Doppler examination is necessary before use of the KITBIFLEX in order to ensure there is no associated severe arterial disease.

• For an ABI between 0.6 and 0.9: the KITBIFLEX must be used with caution, closely monitored, under the supervision of a doctor.

• In the event of discomfort, significant hindrance, pain, variation in limb volume, abnormal sensations or change in colour of the extremities, remove the device and consult a healthcare professional.

• Follow the instructions provided by the health care professional who prescribed or supplied the product, complying with the application method and duration of wear indicated.

• Store at room temperature, preferably in the original packaging.

• Do not apply the product directly to an open wound without a dressing.

• Fully cover any infected or open skin wounds before using the bandage.

• For some conditions (or situations) like venous thrombosis, bandages should be used in combination with anticoagulants; please refer to the advice provided by your health care professional.

This is a medical device supplied on medical prescription. Application requires intervention by a health care professional or a person trained in the application of bandages. Read the instruction leaflet carefully.

Undesirable side-effects

This device may cause a skin reaction such as pressure marks, redness, itching, burning or blisters. Any serious incidents related to the product should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Care/Maintenance

Machine wash, without detergent or bleach.

Bandages No. 1 and No. 2:

Washable 3 times at 40°C, normal wash programme. Dry flat, away from direct heat. Do not bleach. Do not tumble dry. Do not iron. Do not dry clean. Can be autoclaved.

To keep track of the number of washes, tick the boxes on the bandage with a permanent marker suitable for contact with the skin.

Disposal

Dispose of in accordance with local regulations.

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En caso necesario, para reducir el riesgo de deslizamiento, enrollar el enganche largo alrededor de la parte superior de la pantorrilla (ver dibujo ©) sin tensión. El enganche debe rodear todo el miembro y fijarse sobre sí mismo dejando visible el logo Biflex.
THUASNE UK Ltd**
Thuasne SAS – SIREN/RCS Nanterre 542 091 186